BPH CATHETER SYSTEM

Instructions for Use

Table of Contents

1	DEVICE DESCRIPTION	3
1.1	BPH PRODUCT	3
1.2	OPTILUME BPH CATHETER SYSTEM	3
1.2.1	<i>Optilume BPH Catheter</i>	3
2	DRUG COATING	5
3	INTENDED USE	5
4	INDICATIONS FOR USE	5
5	CONTRAINDICATIONS	5
6	WARNINGS	5
7	PRECAUTIONS	6
8	POSSIBLE COMPLICATIONS	8
9	DRUG INFORMATION	8
10	HOW SUPPLIED	9
11	STORAGE	10
12	RECOMMENDED ITEMS	10
13	DIRECTIONS FOR USE	10
13.1	PERI PROCEDURAL MEDICATION	10
13.2	CATHETER PREPARATION	10
13.3	PRE-DILATION	11
13.4	BALLOON RECAPTURE	12
13.5	PRE-DILATION BALLOON COMPLIANCE CHART	12
13.6	OPTILUME BPH DRUG COATED BALLOON DEVICE SIZE SELECTION	12
13.7	DRUG COATED BALLOON CATHETER PREPARATION	13
13.8	DRUG COATED BALLOON DILATION	13
13.9	DRUG COATED BALLOON RECAPTURE	14
13.10	DRUG COATED BALLOON COMPLIANCE CHART	14
14	WARRANTY	15
15	SYMBOLS USED IN THE DEVICE LABELS	16

1 DEVICE DESCRIPTION

The Optilume BPH Catheter System is used to perform a TransUrethral Anterior commissurotomy of the Prostate (TUAP). The distal end of the catheter has a semi-compliant inflatable double lobe balloon that is coated with a proprietary coating containing the active pharmaceutical paclitaxel. The drug coating covers the working length of the balloon body. The Optilume BPH System is comprised of a pre-dilation uncoated balloon catheter and a dilation drug coated balloon (DCB) catheter.

1.1 Optilume BPH Product

The Optilume BPH Catheter device sizes and part numbers are provided in Table 1.

Table 1: Part Number Identification of the Catheter

Part Number	Description
1137-30035B	Optilume BPH Prostatic Pre-dilation Catheter 30x35
1113-30035B	Optilume BPH Prostatic Dilation DCB Catheter 30x35
1113-30040B	Optilume BPH Prostatic Dilation DCB Catheter 30x40
1113-30045B	Optilume BPH Prostatic Dilation DCB Catheter 30x45

1.2 Optilume BPH Catheter System

The Optilume BPH Catheter System consists of the following:

- Optilume BPH Prostatic Pre-dilation catheter (Step 1)
- Optilume BPH Prostatic DCB Catheter (Step 2)

1.2.1 Optilume BPH Catheter

The Optilume BPH Catheter is a single inflation lumen balloon catheter that terminates in an atraumatic tip. The folded balloon has a 14.5Fr profile. The device is inserted through the outer sheath of a rigid cystoscope and can be visualized in a side-by-side fashion with a rigid cystoscope.


Figure 1: Optilume BPH Catheter

Several balloons lengths will be offered to fit different prostate sizes as shown in Table 2.

Table 2: Balloon Size Matrix & Catheter Profile

Balloon Diameter		Balloon Treatment Length			Catheter Profile
mm	Fr	mm			
30	90	35	40	45	14.5 Fr


Figure 2: Inflated Pre-dilation Balloon Catheter with Balloon Protector Sheath Retracted


Figure 3: Inflated Drug Coated Balloon Catheter with Balloon Protector Sheath Retracted

The Optilume BPH Catheter System is used to perform a TransUrethral Anterior commissurotomy of the Prostate (TUAP). It has a double lobe design with a neck separating the two sections. The balloon neck is reinforced preventing diameter growth during the inflation process.

The double lobe design allows the balloon neck to seat in the bladder neck during inflation and helps prevent migration of the balloon into the bladder. The distal lobe of the balloon inflates in the bladder and aids in anchoring the device, while the proximal lobe of the balloon is positioned in the prostatic urethra to dilate the prostate and create the anterior commissurotomy. Both the pre-dilation and the drug coated balloon catheters are identical, the only difference being the DCB is paclitaxel drug coated. The drug coated balloon body and a portion of the proximal cone are coated with the active pharmaceutical ingredient (API) paclitaxel as shown in Figure 3.

The proximal end of the catheter is tapered to allow the attachment of a Tuohy-Borst valve and stopcock for inflation. The balloon is inflated with an inflation device through the Tuohy-Borst connector.

All devices and components are provided sterile and intended for single use only. The catheters are individually packaged in a double pouch configuration.

- The outer pouch is a foil pouch. The interior and exterior of the foil pouch are NON-STERILE and contains a desiccant. The desiccant packs are to be discarded. DO NOT open the foil pouch over the sterile field.
- Inside the foil pouch is a Tyvek pouch. The contents inside the Tyvek pouch are sterile but the external part of the Tyvek pouch is non-sterile.

2 DRUG COATING

The drug coating consists of the active pharmaceutical ingredient paclitaxel and excipients. The drug coating is evenly distributed across the balloon surface as shown in Figure 3 and is the same for all balloon sizes. The key functional characteristic of the drug coating is to allow for release of paclitaxel to the prostatic urothelium during balloon inflation. The table below shows the nominal dose of paclitaxel per balloon size.

Table 3: Nominal Paclitaxel Dose per Balloon [mg]

Balloon Diameter	Balloon Length (mm)		
mm (Fr)	35	40	45
30 (90)	11.4mg	12.5mg	13.7mg

3 INTENDED USE

The Optilume BPH Catheter System is intended to perform a TransUrethral Anterior commissurotomy of the Prostate (TUAP) in adult men with obstructive urinary symptoms associated with Benign Prostatic Hyperplasia (BPH).

4 INDICATIONS FOR USE

The Optilume BPH Catheter system is indicated for the treatment of adult men with obstructive urinary symptoms associated with Benign Prostatic Hyperplasia (BPH).

5 CONTRAINDICATIONS

The Optilume BPH Catheter is contraindicated for use in:
Patients with known hypersensitivity to paclitaxel or structurally related compounds.

6 WARNINGS

- The Optilume BPH Catheter is supplied sterile for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilizing could increase the risk of patient infection and risk of compromised device performance.
- Do not use this device if there is infection in the urinary tract.
- The Optilume BPH catheter should be used only by physicians who are experienced and knowledgeable with the clinical and technical aspects of balloon dilation.
- Prior to use of the Optilume BPH Catheter, physicians should read and understand the instructions for use. Failure to follow the indications, contraindications, restrictions, warnings and precautions may

result in complications.

- Do not use after the “Use By” date.
- Never use air or any gaseous medium to inflate the Optilume BPH Catheter. During use, the Optilume BPH catheter should be manipulated under direct visualization via cystoscopy.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Resistance may cause damage to device or urethra. Carefully withdraw the catheter.
- The BPH DCB contains paclitaxel, a known genotoxin. Men should have protected sex (wear a condom) or abstain for 30 days post treatment.
- Monitor for signs of anaphylaxis or hypersensitivity to paclitaxel.
- Men with sexual partners of childbearing potential should use highly effective contraceptive for at least 6 months post-treatment.
- Adverse reaction to paclitaxel and symptoms observed derived primarily from IV infusion studies of the drug in treating cancer subjects include:
 - Chromosomal abnormalities and the risk of cancer
 - Fetal harm when a pregnant woman is exposed
 - Anaphylaxis and hypersensitivity with paclitaxel
 - Inhibition of the healing of the urethra post procedure
 - Myelosuppression including: neutropenia, leukopenia, thrombocytopenia, anemia
 - Arrhythmia
 - Peripheral neuropathy
 - Myalgia or Arthralgia
 - Alopecia
 - Hypotension
 - Nausea, vomiting or diarrhea
 - Elevated bilirubin, ALP and AST
 - Potential effect on the liver and kidneys is unknown and has not been studied

The amount of the paclitaxel delivered locally during the Optilume BPH procedure is much lower than a single dose of systemic chemotherapy provided to cancer patients. Unlike the chemotherapy application, the drug from the Optilume BPH is deposited in the urethra minimizing systemic exposure. In previous clinical studies using the Optilume BPH DCB, a miniscule amount of paclitaxel was observed systemically for up to 10 hours post-procedure (C_{max} of 0.25 ng/mL immediately post-procedure with the average concentration decreasing below quantifiable levels at 10 hours).

7 PRECAUTIONS

- Carefully inspect the packaging prior to use. Do not use the catheter if it is damaged or if the size, shape or condition is unsuitable for the intended procedure.
- Balloon catheters are intended for use by physicians trained and experienced in techniques for balloon catheter dilation.
- To ensure proper regulation of balloon pressure, use of a balloon catheter inflation device with pressure gauge is recommended.
- Do not immerse or wipe the drug coated catheter with any fluid as the integrity of the drug coating may be damaged or compromised. Replace the drug coated catheter if the balloon has come into contact with fluids prior to use.
- Use dry sterile gloves or dry gauze pads to handle the drug coated catheter prior to use. Care should be taken to minimize contact with the drug coated portion of the device.
- Never inflate the catheter outside the body or prior to reaching the prostatic urethra as it may disrupt the coating integrity.

- Always inflate with a sterile liquid. Never inflate with air, carbon dioxide or any other gas. The balloon catheters should not be inflated in excess of the rated burst pressure (RBP). Do not overinflate the balloon. Inflation in excess of rated burst pressure may cause the balloon to rupture.
- For proper drug delivery to the prostatic urethra, allow the coating to hydrate in situ for approximately 1 minute prior to inflation and maintain inflation of the catheter for a minimum of 5 minutes. To optimize the anterior commissurotomy and drug delivery, longer inflation times > 5 minutes may be performed at the discretion of the operator.
- If the balloon catheter has a failure prior to, or during inflation, replace the balloon catheter and inflate per procedure. If failure of the Optilume drug coated catheter is after inflation to RBP, do not repeat/replace with a second drug coated balloon catheter procedure.
- Aspirate the balloon completely before withdrawing the deflated balloon and removing the device from the urethra. Empty the inflation device and reapply the vacuum to ensure full aspiration of the balloon fluid. Do not use excessive force to withdraw the balloon. Using excessive force to withdraw the balloon can inflict trauma to tissue.
- Handle and dispose of the used device in accordance with accepted medical practice and applicable local regulations for biohazard waste.
- During the healing phase, the use of a Foley catheter may be required for several days. If a post procedural catheter is used, traction should be applied to firmly seat the Foley balloon against the bladder neck. If moderate bleeding is observed through the catheter drainage line, it is recommended to flush the bladder with 1 liter of sterile saline prior to discharging the patient.
- Safety and effectiveness data has not been established during the clinical study to support the treatment of BPH in patients with:
 - Penile implant or stent
 - Prior minimally invasive intervention (TUNA, Balloon, Microwave, Rezum, Urolift) or surgical intervention of the prostate.
 - Prostate Cancer
 - Previous Radical prostatectomy
 - Pelvic radiation
 - Visible hematuria
 - Neurogenic bladder or sphincter abnormalities or neurological disorders that might affect bladder or sphincter function
 - Previous or current diagnosis of urethral strictures, bladder neck contracture or detrusor muscle spasms
 - Patients that have used beta blockers, antihistamines, anticonvulsants, or antispasmodics within 1 week prior to treatment unless there is evidence of stable dosing for last 6 months
 - Patients that have used alpha blockers, antidepressants, anticholinergics, androgens, daily tadalafil or gonadotropin-releasing hormonal analogs within 3 weeks prior to treatment
 - Patient that have used 5-alpha reductase inhibitor within 6 months prior to treatment
 - Incidence of spontaneous urinary retention within 6 months prior to baseline assessment
 - Post-void residual volume > 250 ml or catheter dependent bladder drainage
 - Overactive bladder (OAB) or urge incontinence
 - Known poor detrusor muscle function (e.g. Qmax < 5 ml/sec)
 - Current bladder stones or prostatic calculi
 - Biopsy of prostate within 30 days prior to procedure or planned within 30 days following the procedure
 - History of cancer in non-genitourinary system which is not considered cured (except basal cell or squamous cell carcinoma of the skin).
 - History of clinically significant comorbidities or presence of unstable conditions (e.g.

cardiovascular, lung, renal [serum creatinine > 2.0 mg/dl], hepatic, bleeding disorders, or metabolic impairment)

- Expected life expectancy < one year
- Current treatment with anti-coagulants (e.g., warfarin or enoxaparin) or anti-platelet medications other than aspirin (e.g., clopidogrel)
- Intravesical prostatic protrusion (IPP) > 1 cm
- Current uncontrolled diabetes (hemoglobin A1c > 7%)

8 POSSIBLE COMPLICATIONS

Possible complications associated with the use of the Optilume BPH Catheter System may include, but are not limited to:

- Pain and tenderness
- Bladder spasm from Foley catheter placement
- Tissue Trauma in surrounding structures including urethral damage
- Hematuria
- Drug reactions, allergic reaction to contrast medium used during diagnostic urethrogram
- Urinary Tract Infection (UTI)
- Tissue perforation
- Incontinence
- Dysuria
- Fever
- Urinary retention

9 DRUG INFORMATION

○ MECHANISM OF ACTION

The Optilume BPH drug coating contains paclitaxel, an anti-mitotic pharmaceutical agent that specifically binds to and stabilizes microtubules. Paclitaxel has been reported to inhibit smooth muscle cell and fibroblast proliferation and migration as well as secretion of extracellular matrix. The combination of these effects may result in the inhibition of prostatic hyperplasia.

○ DRUG INTERACTIONS

Formal drug interaction studies have not been conducted for the Optilume BPH catheter. The respective instructions for use for all drugs used in conjunction with the Optilume BPH catheter should be consulted for interactions with paclitaxel.

Consideration should be given to the potential for systemic and local drug interactions in the prostate in a patient who is taking a drug with known interactions to paclitaxel or when deciding to initiate drug therapy in a patient who has been treated with the BPH Prostatic DCB.

The metabolism of paclitaxel is catalyzed by cytochrome P450 isoenzymes CYP2C8 and CYP3A4, and it is a substrate of P-glycoprotein. Potential drug interactions may occur with any drug that affects these isoenzymes. In the absence of formal drug interaction studies, caution should be exercised when administering paclitaxel.

○ CARCINOGENICITY, GENOTOXICITY AND REPRODUCTIVE TOXICOLOGY

No long-term studies have been performed to evaluate the carcinogenic potential of the drug paclitaxel or of the Optilume BPH Prostatic DCB, and there are no adequate and well-controlled studies published in pregnant women or in men intending to father children. Paclitaxel inhibits cell proliferation by interacting with microtubules, and one consequence is the loss of whole chromosomes during cell division. This indirect action is consistent with positive responses in vitro and in vivo micronucleus genotoxicity assays, which detect DNA fragments. Positive results have also been reported for chromosomal aberrations in primary human lymphocytes. It is not known whether paclitaxel has a separate direct action on DNA in the generation of DNA strand breaks or fragments. It is negative in assays for gene mutation, including salmonella and CHO/HPRT.

Studies performed in rats and rabbits receiving IV paclitaxel during organogenesis revealed evidence of maternal toxicity, embryotoxicity, and fetotoxicity at dosages of 1 and 3 mg/kg, respectively [approximately 5 and 15 times the dose provided by the Optilume BPH DCB coated with 13.7 mg paclitaxel (30mm x 45mm balloon) adjusted for body weight]. No teratogenicity was observed in gravid rats receiving daily IV paclitaxel doses of 1 mg/kg [a daily dose of approximately 5 times the dose of the Optilume BPH DCB (30mm x 45mm), adjusted for bodyweight].

The treating physician should balance the potential medical benefits of the Optilume BPH Catheter against these genotoxic and reproductive risks. **WARNING:** The Optilume BPH Prostatic DCB catheter contains paclitaxel, a known genotoxin. Men should have protected sex (wear a condom) for 30 days post treatment.

○ CARDIOVASCULAR METADATA ANALYSIS

A recent analysis of study data suggests the possibility of an increased risk of death resulting from the use of paclitaxel-coated vascular devices.

Combined data from several different clinical studies appeared to show that patients with peripheral arterial (vascular) disease who were treated with paclitaxel-coated balloons and paclitaxel-releasing stents were at a higher risk of death. This higher risk of death was observed at time points at least two years after treatment with the paclitaxel-containing devices. The reason for this increased risk of death or how strong the risk might be is unknown. Whether use of these devices in future patients with peripheral arterial disease will demonstrate similar additional risk of death is also unknown.

We also do not know whether the risk of death observed in patients treated with paclitaxel-coated devices inside their blood vessels means there is also a risk of death for patients treated in other areas of the body, such as the prostate. Patients receiving the Optilume BPH DCB will be treated with a different paclitaxel-coated balloon for a different condition (BPH) in a different part of the body (the prostate). Nevertheless, based on the known information, one of the potential risks of this treatment is an increased risk of death due to the exposure to paclitaxel during the treatment.

However, unlike the cardiovascular application, the drug is deposited in the prostate and not in the blood, although a small amount of drug can diffuse through the prostate into the blood. Previous urological studies measuring the drug in the blood with the Optilume BPH DCB showed that the majority of the drug was removed in the urine and very little drug got into the blood stream (systemic exposure).

10 HOW SUPPLIED

Each Optilume BPH Catheter is supplied STERILE for single use only (ethylene oxide sterilization). Each

balloon catheter is in a double pouch packaging system (foil and Tyvek pouches) contained within a single unit box. The foil and external surface of the Tyvek pouch is nonsterile. Only the contents within the Tyvek pouch are sterile.

11 STORAGE

The Optilume BPH Catheter should be stored at room temperature between 15°C - 30°C in a dry location in its original packaging. The device should be used prior to the “Use by” date on the packaging.

12 RECOMMENDED ITEMS

Prepare the following items using sterile technique:

- Lubricious surgical gel to aid with device insertion
- 60 ml inflation device with pressure gauge
- Three-Way Stopcock (provided with inflation device)
- Tuohy-Borst adapter
- Rigid cystoscope (minimum sheath size 19.5F)
- Sterile saline
- Sterile contrast media (if using fluoroscopy)

13 DIRECTIONS FOR USE

13.1 Peri Procedural Medication

It is recommended that physicians follow guidelines for pre-procedure medications and preparation for a cystoscopic urethral procedure, including the administration of a pre-procedure antibiotic as appropriate. Oral NSAIDs are also recommended to be given prior to the procedure.

If a urinary tract infection (UTI) is present at the time of treatment, the patient must be treated until the infection is cleared before the Optilume dilation procedure can take place.

13.2 Catheter Preparation

Predilation of the prostatic urethra should be completed with the Optilume BPH Pre-dilation Catheter (1137-30035) to create an anterior commissurotomy prior to treatment with the Optilume BPH drug coated balloon catheter.

Evacuate air from the predilation balloon catheter. The balloon lumen of the catheter contains air and the air must be displaced to make certain that only liquid fills the balloon while the catheter is in the urethra.

1. Open the sterile inflation device package and remove the 3-way stopcock and attach it to the Tuohy-Borst adapter.
2. Connect the Tuohy-Borst adapter to the proximal catheter shaft. Ensure that the catheter shaft does not interfere with the 3-way stopcock operation.
3. Fill the inflation device half-way with sterile saline and attach it to the stopcock. Turn the stopcock so fluid can flow between the inflation device and the balloon catheter.
4. With the inflation device pointing downwards, draw back plunger to full volume of syringe (this creates maximum negative pressure and allow air to evacuate above the fluid level) and hold until no air bubbles can be seen coming out of the saline in the syringe. Repeat as needed to purge the air from the catheter.

5. With catheter preparation complete, disconnect the Tuohy-Borst adapter from the catheter shaft.
6. Fill the inflation device with 50cc of normal saline and purge air from the line.

13.3 Pre-dilation

1. Select and prepare the Optilume BPH Predilation Prostatic Balloon catheter per Table 4 and Section 13.2.
2. Assemble and advance the rigid cystoscope (minimum sheath size 19.5Fr) through the urethra and into the bladder. Remove the optics and bridge (if applicable), leaving an open pathway through the rigid sheath into the bladder.
3. Remove the Balloon Protector Sheath (BPS) from the balloon catheter
4. Insert the balloon catheter through the rigid cystoscope sheath and into the bladder.
5. Slide the outer rigid cystoscope sheath out of the patient and over the balloon catheter shaft, while maintaining the balloon catheter in position.
6. Reassemble the rigid cystoscope sheath with the bridge and optics.
7. Attach the Tuohy-Borst adapter with inflation device to the proximal catheter shaft.
8. Insert the reassembled cystoscope transurethrally up to the external sphincter. The Optilume BPH device is used side-by-side with a rigid cystoscope.
9. Locate the external sphincter with the cystoscope and position the tip of the cystoscope so visualization of the external sphincter can be maintained throughout the procedure.
10. Adjust the position of the balloon by pushing/pulling the catheter shaft until the proximal balloon is located in the prostatic urethra. Proper positioning is achieved when the proximal balloon bond (where the balloon attaches to the catheter shaft) is visualized with the cystoscope and positioned within the external sphincter per Figure 4.


Figure 4: Positioning of the Balloon Within the Prostatic Urethra

11. With the balloon properly positioned, predilate the prostate using the inflation device. Inflate slowly, approximately 1 ATM every 30 seconds while maintaining traction to seat the balloon at the bladder neck and prevent proximal migration into the bladder. During inflation monitor the location of the proximal balloon bond with the cystoscope. If the proximal balloon bond is no longer visible within the external sphincter or you believe the treatment balloon has slipped into the bladder, deflate and reposition the balloon as above and repeat the dilation process.

NOTE: A blue marker is placed approximately 8mm proximal to the proximal balloon bond on the catheter shaft. This marker can be utilized to aid in balloon positioning and determining if the balloon has slipped into the bladder.

- Once the balloon is inflated check for migration into the bladder by grabbing the catheter shaft and pressing it gently towards the bladder. If the catheter shaft moves freely, the balloon has likely migrated into the bladder. If migration is observed, deflate the balloon and reposition the catheter as above and repeat the dilation process. Do not repeat this procedure more than 3 times

Warning: If you observe migration of the balloon into the external sphincter, stop and reposition.

- Continue slow inflation until the anterior prostatic commissure is separated indicated by a sudden drop in pressure on the inflation device. Do not exceed the RBP (4 atm) of the balloon.
- After dilation, deflate the balloon and visually confirm the anterior commissurotomy by advancing the cystoscope into the prostatic urethra.
- Do not exceed 3 inflation cycles to rated burst pressure. If a commissurotomy is not achieved after 3 inflation attempts to RBP, recapture (see section 13.4 below) the Pre-dilation balloon and proceed to dilation with the DCB.

13.4 Balloon Recapture

- With the commissurotomy achieved, remove the cystoscope. Note: the Pre-dilation balloon initiates the anterior commissurotomy while the DCB dilation completes the commissurotomy.
- Deflate balloon by applying aspiration with an inflation device. Once the inflation device is filled with fluid, detach from the catheter, purge fluid into a reservoir, and reattach to the catheter. The inflation device needs to be emptied and vacuum reapplied to ensure full aspiration of the balloon fluid.
- When the balloon is completely deflated, maintain vacuum and gently pull on the catheter shaft to withdraw the catheter from the patient's body.

13.5 Pre-dilation Balloon Compliance Chart

30mm(90Fr) x 35mm

(ATM) Pressure	kPa		(mm) Balloon
1.0	101		28.70
2.0	203	Nominal	30.15 (90Fr)
3.0	304		31.36
4.0	405	RBP	32.33

13.6 Optilume BPH Drug Coated Balloon Device Size Selection

Use a transrectal ultrasound (TRUS) to determine the prostatic urethral length (PUL), volume, and height of the prostate. Select the corresponding balloon length for the treatment from the table below.

Table 4: BPH Drug Coated Balloon Selection Matrix

Prostate Volume & Height	Prostatic Urethra Length		
	PUL 37-42mm	PUL 42-47mm	PUL >47mm
Vol 20 to 80g Prostate Height ≥ 20mm	30x35mm	30x40mm	30x45mm

- A smaller balloon size can always be used at the physician's discretion.

13.7 Drug Coated Balloon Catheter Preparation

Repeat steps listed in section 13.2.

13.8 Drug Coated Balloon Dilation

1. Select and prepare the Optilume BPH Drug Coated Balloon catheter per Table 4 and Section 13.2.
2. Assemble and advance the rigid cystoscope (minimum sheath size 19.5Fr) through the urethra and into the bladder. Remove the optics and bridge (if applicable), leaving an open pathway through the rigid sheath into the bladder.
3. Remove the Balloon Protector Sheath (BPS) from the balloon catheter.
4. Insert the balloon catheter through the rigid cystoscope and into the bladder.
5. Slide the outer rigid cystoscope sheath out of the patient and over the balloon catheter shaft while maintaining the balloon catheter in position.
6. Reassemble the rigid cystoscope sheath with the bridge and optics.
7. Attach the Tuohy-Borst adapter with inflation device to the proximal catheter shaft.
8. Insert the reassembled cystoscope transurethraly up to the external sphincter. The Optilume BPH device is used side-by-side with a rigid cystoscope.
9. Locate the external sphincter with the cystoscope and position the tip of the cystoscope so visualization of the external sphincter can be maintained through the procedure.
10. Adjust the position of the balloon by pushing/pulling the catheter shaft until the proximal balloon bond (where the balloon attaches to the catheter shaft) is visualized with the cystoscope and positioned within the external sphincter per Figure 5. This step should be completed in approximately 1 minute while the coating is hydrating, such that the balloon can be inflated after hydration is complete.


Figure 5: Positioning of the Balloon Within the Prostatic Urethra

11. With the balloon properly positioned, dilate the prostate using the inflation device. Do not exceed the RBP of the balloon.

Warning: If you observe migration of the balloon into the external sphincter, stop, deflate and reposition.

12. Maintain inflation for a minimum of 5 minutes to assure complete anterior commissurotomy.

13.9 Drug Coated Balloon Recapture

With dilation complete, remove the cystoscope.

1. Deflate balloon by applying aspiration with an inflation device. Once the inflation device is filled with fluid, detach from the catheter, purge fluid into a reservoir, and reattach to the catheter. The inflation device needs to be emptied and vacuum reapplied to ensure full aspiration of the balloon fluid.
2. When the balloon is completely deflated, maintain vacuum and gently pull on the catheter shaft to withdraw the catheter from the patient's body.
3. Dispose of properly in accordance with accepted medical practice and applicable local regulations for biohazard waste.

13.10 Drug Coated Balloon Compliance Chart

30mm(90Fr) x 35mm

(ATM) Pressure	kPa		(mm) Balloon
1.0	101		28.70
2.0	203	Nominal	30.15 (90Fr)
3.0	304		31.36
4.0	405	RBP	32.33

30mm(90Fr) x 40mm

(ATM) Pressure	kPa		(mm) Balloon
1.0	101		28.39
2.0	203	Nominal	30.07 (90Fr)
3.0	304		31.49
4.0	405	RBP	32.43

30mm(90Fr) x 45mm

(ATM) Pressure	kPa		(mm) Balloon
1.0	101		28.97
2.0	203	Nominal	30.75 (90Fr)
3.0	304		32.14
4.0	405	RBP	33.07

Caution: The rated burst pressure should not be exceeded. Refer to product label for rated burst pressures. Inflation beyond the rated burst pressure may cause the balloon to rupture. If loss of pressure within the balloon occurs during inflation or if balloon ruptures during dilation, immediately discontinue the procedure. Deflate the balloon carefully and remove from urethra. Do not re-inflate.

14 WARRANTY

Urotronic warrants that reasonable care has been used in the design and manufacture of this product. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties for a particular purpose. Handling, storage, cleaning and sterilization of this device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Urotronic's control directly affect the device and the results obtained from its use. Urotronic's obligation under this warranty is limited to the repair or replacement of this device and Urotronic shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this device. Urotronic assumes no liability with respect to devices reused, reprocessed or re-sterilized and makes no warranties, express or implied, including but not limited for a particular purpose, with respect to such devices.

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15 SYMBOLS USED IN THE DEVICE LABELS

	Quantity of 1 per box
	Caution: Federal law restricts this device to sale by or on the order of a physician
	Indicates the date when the medical device was manufactured
	Do not re-sterilize
	Do not re-use
	Do not use if package is damaged
	Fragile
	Use-by date
	Keep away from sunlight
	Keep Dry
	Manufacturer
	Does not contain latex
	Temperature limit 15°C - 30°C
	Caution: Consult instructions for use
	Sterilized using ethylene oxide

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MDSS GmbH